

# Atropos™ SC

## 0.035" system

Combines strength and optimal control in routine and challenging cases

Designed to treat fibrotic and calcified lesions

Broad size range:  
Diameter: 3.0 - 10.0mm  
Length: 20 - 200mm

## Semi-Compliant OTW

Engineered for AV access

### Excellent pushability


Seamless tip-to-shaft design for smooth push force transmission

### Lowest tip profile 0.038"

Facilitates lesion crossing with ease

### Slider™ hydrophilic coating

Latest hydrophilic coating technology to facilitate easy maneuvering through complex anatomy and lesions



# Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	20	30	40	60	80	100	120	150	200				
3.0	V	V	V	V	-	-	-	-	-	40	5	6	14
4.0	V	V	V	V	-	-	-	-	-	40	5	6	14
5.0	V	V	V	V	-	-	-	-	-	40	5	6	14
6.0	V	V	V	V	-	-	-	-	-	40	6	6	14
7.0	V	V	V	V	-	-	-	-	-	40	6	6	14
8.0	V	V	V	V	-	-	-	-	-	40	6	6	14
3.0	V	V	V	V	V	V	V	V	V	70	5	6	14
4.0	V	V	V	V	V	V	V	V	V	70	5	6	14
5.0	V	V	V	V	V	V	V	V	V	70	5	6	14
6.0	V	V	V	V	V	V	V	V	V	70	6	6	14
7.0	V	V	V	V	V	V	V	V	V	70	6	6	14
8.0	V	V	V	V	V	V	V	V	—	70	6	6	14
9.0	V	V	V	V	V	V	V	—	—	70	7	6	14
10.0	V	V	V	V	V	—	—	—	—	70	7	6	14
3.0	V	V	V	V	V	V	V	V	V	90	5	6	14
4.0	V	V	V	V	V	V	V	V	V	90	5	6	14
5.0	V	V	V	V	V	V	V	V	V	90	5	6	14
6.0	V	V	V	V	V	V	V	V	V	90	6	6	14
7.0	V	V	V	V	V	V	V	V	V	90	6	6	14
8.0	V	V	V	V	V	V	V	V	—	90	6	6	14
9.0	V	V	V	V	V	V	V	—	—	90	7	6	14
10.0	V	V	V	V	V	—	—	—	—	90	7	6	14
3.0	V	V	V	V	V	V	V	V	V	150	5	6	14
4.0	V	V	V	V	V	V	V	V	V	150	5	6	14
5.0	V	V	V	V	V	V	V	V	V	150	5	6	14
6.0	V	V	V	V	V	V	V	V	V	150	6	6	14
7.0	V	V	V	V	V	V	V	V	V	150	6	6	14
8.0	V	V	V	V	V	V	V	V	—	150	6	6	14
9.0	V	V	V	V	V	V	V	—	—	150	7	6	14
10.0	V	V	V	V	V	—	—	—	—	150	7	6	14



## BrosMed Medical Co., Ltd.

Add: 15th building, SMEs Park, Songshan Lake Hi-Tech Zone, 523808, Dongguan, Guangdong, China  
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016  
 E-mail: sales@brosmed.com www.brosmed.com

## BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard The Netherlands  
 Tel: +31 (0) 850 140 901  
 E-mail: cs@brosmed.com  
 Web: www.brosmed.com

US FDA K160256



2797  
 GRA-6022/Rev.01