

Balloon Inflation Devices



BM-BID- I



BM-BID- II



BM-BID- III

INTRODUCTION

This product is used to control the inflation and deflation of various balloon catheters in minimally invasive interventional therapy or diagnostic surgery.

PRODUCT FEATURES

- Ergonomic grip provides comfort and convenience
- Accurate pressure control
- Clock face is clear and easy to read

ORDERING INFORMATION

Catalogue No.	Volume (ml)	Pressure (ATM)
BM-BID-I30	20	30
BM-BID-I30L	30	30
BM-BID-II26	20	26
BM-BID-II30	20	30
BM-BID-III30	20	30