

Castor NC

0.014" system

Excellent performance in routine and challenging cases

Designed to treat resistant lesions


Optimized balloon strength
RBP up to 22atm

Non-Compliant OTW

Designed to treat BTK lesions

Optimized crossing profile

Best in class crossing capability for diffused and challenging lesions



Specifications

Balloon Diameter (mm)	Balloon Length (mm)									Catheter Working Length (cm)	Sheath Compatibility (F)	Nominal Pressure (atm)	Rated Burst Pressure (atm)
	10	15	20	30	40	60	80	100	150				
2.0	V	V	V	V	V	V	V	V	V	70	4	12	22
2.5	V	V	V	V	V	V	V	V	V	70	4	12	22
3.0	V	V	V	V	V	V	V	V	V	70	4	12	22
3.5	V	V	V	V	V	V	V	V	V	70	4	12	22
4.0	V	V	V	V	V	V	V	V	V	70	4	12	22
4.5	V	V	V	V	V	V	V	V	V	70	5	12	20
5.0	V	V	V	V	V	V	V	V	V	70	5	12	20
5.5	V	V	V	V	V	V	V	V	V	70	5	12	20
6.0	V	V	V	V	V	V	V	V	V	70	5	12	20
7.0	V	V	V	V	V	V	V	V	V	70	6	12	18
8.0	V	V	V	V	V	V	V	V	V	70	6	12	18
2.0	V	V	V	V	V	V	V	V	V	90	4	12	22
2.5	V	V	V	V	V	V	V	V	V	90	4	12	22
3.0	V	V	V	V	V	V	V	V	V	90	4	12	22
3.5	V	V	V	V	V	V	V	V	V	90	4	12	22
4.0	V	V	V	V	V	V	V	V	V	90	4	12	22
4.5	V	V	V	V	V	V	V	V	V	90	5	12	20
5.0	V	V	V	V	V	V	V	V	V	90	5	12	20
5.5	V	V	V	V	V	V	V	V	V	90	5	12	20
6.0	V	V	V	V	V	V	V	V	V	90	5	12	20
7.0	V	V	V	V	V	V	V	V	V	90	6	12	18
8.0	V	V	V	V	V	V	V	V	V	90	6	12	18
2.0	V	V	V	V	V	V	V	V	V	150	4	12	22
2.5	V	V	V	V	V	V	V	V	V	150	4	12	22
3.0	V	V	V	V	V	V	V	V	V	150	4	12	22
3.5	V	V	V	V	V	V	V	V	V	150	4	12	22
4.0	V	V	V	V	V	V	V	V	V	150	4	12	22
4.5	V	V	V	V	V	V	V	V	V	150	5	12	20
5.0	V	V	V	V	V	V	V	V	V	150	5	12	20
5.5	V	V	V	V	V	V	V	V	V	150	5	12	20
6.0	V	V	V	V	V	V	V	V	V	150	5	12	20
7.0	V	V	V	V	V	V	V	V	V	150	6	12	18
8.0	V	V	V	V	V	V	V	V	V	150	6	12	18



BrosMed Medical Co., Ltd.

Add: 15th Building, SMEs Venture Park SongShan Lake
 Hi-Tech Industrial Development Zone Dongguan
 523808, GuangDong, China
 Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
 E-mail: sales@brosmed.com www.brosmed.com

BrosMed Medical B.V. (International Sales Office)

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
 Tel: +31 (0) 850 140 901
 E-mail: cs@brosmed.com
 Web: www.brosmed.com

US FDA 510K cleared



2797
 GRA-6025/Rev.01