


Apollo^{NC}

Stronger • Better • Safer

Non-Compliant, High Pressure
Coronary Balloon Catheter



**High Pressure •
Non-Compliant •
6F Guiding Kissing •
3.0% Axial Growth •
0.55% Compliance •**

Apollo is a non-compliant high pressure post-dilatation catheter offering numerous features for enhanced and robust performance :

- Exceptionally minimal axial growth during inflation
- 5F/6F guiding catheter compatible
- Versatile, Optimal balloon diameters and lengths to meet clinical demands (2.0-5.0 mm diameters) (6-30 mm lengths)
- Average balloon burst pressure is 30 atm

Brosmed
embrace life

Balloon Specifications

Balloon Diameter (mm)	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0
Balloon Length (mm)	6, 8, 10, 12, 15, 18, 20, 22, 25, 30
Nominal Pressure	12 atm
Rated Burst Pressure	22 atm (φ2.0-4.0mm) 20 atm (φ4.5-5.0mm)
Average Burst Pressure	30 atm
Compliance	Non-compliant
Balloon Pleats	3
Tip Length (mm)	2.0/2.5
Distal Shaft OD	2.55F/2.6F
Proximal Shaft OD	2.0F
Guiding Catheter Compatibility	5F/6F Compatible, 6F Compatible for kissing
Working Length (cm)	140
Coating	Hydrophilic external surface coating from tip to adjacent guidewire port; Silicone Coating inside guidewire lumen

Ordering Information

Balloon Diameter	Balloon Length / Product Code									
	6mm	8mm	10mm	12mm	15mm	18mm	20mm	22mm	25mm	30mm
ø2.0mm	802-2006	802-2008	802-2010	802-2012	802-2015	802-2018	802-2020	802-2022	802-2025	802-2030
ø2.25mm	802-2206	802-2208	802-2210	802-2212	802-2215	802-2218	802-2220	802-2222	802-2225	802-2230
ø2.5mm	802-2506	802-2508	802-2510	802-2512	802-2515	802-2518	802-2520	802-2522	802-2525	802-2530
ø2.75mm	802-2706	802-2708	802-2710	802-2712	802-2715	802-2718	802-2720	802-2722	802-2725	802-2730
ø3.0mm	802-3006	802-3008	802-3010	802-3012	802-3015	802-3018	802-3020	802-3022	802-3025	802-3030
ø3.25mm	802-3206	802-3208	802-3210	802-3212	802-3215	802-3218	802-3220	802-3222	802-3225	802-3230
ø3.5mm	802-3506	802-3508	802-3510	802-3512	802-3515	802-3518	802-3520	802-3522	802-3525	802-3530
ø3.75mm	802-3706	802-3708	802-3710	802-3712	802-3715	802-3718	802-3720	802-3722	802-3725	802-3730
ø4.0mm	802-4006	802-4008	802-4010	802-4012	802-4015	802-4018	802-4020	802-4022	802-4025	802-4030
ø4.5mm	802-4506	802-4508	802-4510	802-4512	802-4515	802-4518	802-4520	802-4522	802-4525	802-4530
ø5.0mm	802-5006	802-5008	802-5010	802-5012	802-5015	802-5018	802-5020	802-5022	802-5025	802-5030



BrosMed Medical Co., Ltd.

15th Building, SMES Venture Park, Songshan Lake Hi-Tech Industrial
Development Zone, Dongguan523808, Guangdong, China
Tel: +86 (769) 2289 2018 Fax: +86 (769) 2289 2016
E-mail: sales@brosmed.com

BrosMed Medical B.V.

Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Office: +31(0) 850 140 901
E-mail: cs@brosmed.com

510(K) Cleared

CE 2797

www.brosmed.com

GRA-6002(CE)/Rev.10