


Introducer Sheath Sets



PRODUCT FEATURES

- The sheath offers excellent kink resistance and flexibility with smooth transitions between cannula, dilator and guide wire
- Smooth tapered tip results in low insertion force, avoids wrinkling and reduces puncture site trauma
- Advanced engineering and special design of hemostatic valve maintains good hemostasis and facilitates smooth operation for guide wire and catheter
- Radiopaque sheath material facilitates visibility under fluoroscopy
- Various specifications and kit configurations available

Exclusively for Radial introducer Sheath Set

- 20G I.V. cannula in set to reduce trauma, minimize blood loss and to provide fast hemostasis with a tourniquet
- Smooth kink-resistant guide wire for convenient handling

KIT CONTENTS

- Introducer sheath, dilator, introducer needle, kink-resistant guide wire, scalpel and syringe

DESIGNED FOR RADIAL ACCESS

Catalogue No.	Size (Fr.)	Length (cm)	Polymer Jacket Guidewire	I.V. Cannula
BM-RA0407	4F	7	0.025" × 45cm Straight	20G × 32mm
BM-RA0507	5F	7		
BM-RA0607	6F	7		
BM-RA0707	7F	7		
BM-RA0411	4F	11		
BM-RA0511	5F	11		
BM-RA0611	6F	11		
BM-RA0711	7F	11		
BM-RA0416	4F	16		
BM-RA0516	5F	16		
BM-RA0616	6F	16		
BM-RA0716	7F	16		

DESIGNED FOR RADIAL ACCESS

Catalogue No.	Size (Fr.)	Length (cm)	Nitinol guide wire	Introducer Needle
BM-RA0407B	4F	7	0.018" × 45cm Straight	21G × 40mm
BM-RA0507B	5F	7		
BM-RA0607B	6F	7		
BM-RA0707B	7F	7		
BM-RA0411B	4F	11		
BM-RA0511B	5F	11		
BM-RA0611B	6F	11		
BM-RA0711B	7F	11		

DESIGNED FOR FEMORAL ACCESS

Catalogue No.	Size (Fr.)	Length (cm)	Nitinol guide wire	Introducer Needle
BM-FA0409A	4F	9	0.021" × 45cm Angled	20G × 38mm
BM-FA0511A	5F	11	0.035" × 45cm Angled	18G × 65mm
BM-FA0611A	6F	11		
BM-FA0711A	7F	11		
BM-FA0811A	8F	11		
BM-FA0911A	9F	11	0.035" × 60cm Angled	
BM-FA0419A	4F	19		
BM-FA0523A	5F	23		
BM-FA0623A	6F	23		
BM-FA0723A	7F	23		
BM-FA0823A	8F	23	0.038" × 45cm Angled	
BM-FA0923A	9F	23		
BM-FA0511D	5F	11		
BM-FA0611D	6F	11		
BM-FA0711D	7F	11		
BM-FA0811D	8F	11	0.038" × 60cm Angled	
BM-FA0911D	9F	11		
BM-FA0419D	4F	19		
BM-FA0523D	5F	23		
BM-FA0623D	6F	23		
BM-FA0723D	7F	23	0.038" × 60cm Angled	
BM-FA0823D	8F	23		
BM-FA0923D	9F	23		

Catalogue No.	Size (Fr.)	Length (cm)	Stainless steel guide wire	Introducer Needle
BM-FA0409	4F	9	0.021" × 45cm Angled	20G × 38mm
BM-FA0511	5F	11	0.035" × 45cm Angled	18G × 65mm
BM-FA0611	6F	11		
BM-FA0711	7F	11		
BM-FA0811	8F	11		
BM-FA0911	9F	11	0.035" × 60cm Angled	
BM-FA0419	4F	19		
BM-FA0523	5F	23		
BM-FA0623	6F	23		
BM-FA0723	7F	23		
BM-FA0823	8F	23	0.038" × 45cm Angled	
BM-FA0923	9F	23		
BM-FA0511B	5F	11		
BM-FA0611B	6F	11		
BM-FA0711B	7F	11		
BM-FA0811B	8F	11	0.038" × 60cm Angled	
BM-FA0911B	9F	11		
BM-FA0419B	4F	19		
BM-FA0523B	5F	23		
BM-FA0623B	6F	23		
BM-FA0723B	7F	23	0.038" × 60cm Angled	
BM-FA0823B	8F	23		
BM-FA0923B	9F	23		